



**KEYSURGICAL®**



8101 Wallace Road  
Eden Prairie, MN 55344 USA

Tel. 800.541.7995 952.914.9789  
Fax 952.914.9866

keysurgical.com info@keysurgical.com

## **KEY SURGICAL DISINFECTANTS AND CLEANERS**

Catalog pages 100 - 109

## DISINFECTANTS AND CLEANERS

Key Surgical Disinfectants and Cleaners are designed to enable the quick, easy and efficient disinfection of surfaces and surgical equipment.

Sani-Cloth® germicidal wipes offer a highly effective line of defense against microorganisms. The durable non-woven disposable cloths are designed to clean and disinfect instruments and surfaces such as stainless steel, glass, tables, carts, counters, cabinets, telephones, and other hard nonporous surfaces. Available as 2, 5, or 10 minute disinfectants. All are intended to be rotated to ensure that microorganisms don't become immune to any one formula.

Our line of Cavicide® disinfectant cleaner is ideal for use on all nonporous surfaces including plastic, stainless steel, glass, tile, and more. This unique non-caustic formula offers a 3 minute kill time for TB and a 2 minute kill time for HBV, HCV, HIV-1, MRSA, and VRE. Highly effective at removing mold and mildew. Powerful enough to kill on contact yet gentle enough to be used on ultrasonic instruments and as a decontaminant for other sensitive medical instruments, all while being fully biodegradable. The Cavicide solution is ready to use and requires no rinsing. It is formulated for the control of cross contamination and can be used in 3 different ways: an immersion solution, a surface disinfectant cleaner, or as an ultrasonic solution.

From spray on foams to concentrated multi-purpose detergents, our line of Enzymatic Detergents pre-clean, prevent corrosion, reduce the need for scrubbing and extend the life of your instruments. Safe on metals, plastics and rubbers, our detergents are non-hazardous, non-corrosive and biodegradable.

For a deeper level of cleaning that's still gentle on instruments, EfferZyme Cleaning Tablets help remove blood and other soils that may be resistant to chemical detergents. Wypalls™ are soft, absorbent and excellent for use as general purpose cleaners. Our Polyester Wipers, are made from 100% medium-weight polyester and are perfect for cleanroom use.

Sani-Cloth® Surface Wipes

CaviCide Disinfectant Cleaner

Enzymatic Detergent

EfferZyme Cleaning Tablets

Multi-Purpose Detergent

Kimberly-Clark Wypall™

Polyester Wiper



# Sterile Processing Products – Disinfectants and Cleaners



**SUPER SANI-CLOTH®\***

A nonwoven disposable cloth for use on hard, nonporous surfaces. Cleans, disinfects and deodorizes. 55% alcohol content. Effective in 2 minutes or less against 22 microorganisms. Latex-free.

\* Super Sani-Cloth, Sani-Cloth Plus and Sani-Cloth HB are intended to be used in rotation to ensure that microorganisms don't become resistant to any one type of wipe.

Canister of Extra Large Germicidal Surface Wipes. 65 wipes/can, 8" x 14" wipes.

| Number    | Quantity |
|-----------|----------|
| DI-1065   | 1/pkg    |
| DI-1065-6 | 6/case   |

Canister of Large Germicidal Surface Wipes. 160 wipes/can, 6" x 6.75" wipes.

| Number     | Quantity |
|------------|----------|
| DI-1160    | 1/pkg    |
| DI-1160-12 | 12/case  |

**SANI-CLOTH® PLUS\***

A nonwoven disposable cloth for use on hard nonporous surfaces. Cleans, disinfects and deodorizes. 14.85% alcohol content. Effective in 5 minutes or less against 13 microorganisms. Latex-free.

Canister of Extra Large Germicidal Surface Wipes. 65 wipes/can, 8" x 14" wipes.

| Number    | Quantity |
|-----------|----------|
| DI-2065   | 1/pkg    |
| DI-2065-6 | 6/case   |

Canister of Large Germicidal Surface Wipes. 160 wipes/can, 6" x 6.75" wipes.

| Number     | Quantity |
|------------|----------|
| DI-2160    | 1/pkg    |
| DI-2160-12 | 12/case  |

**SANI-CLOTH® HB\***

A nonwoven disposable cloth for use on hard, nonporous surfaces. Cleans, disinfects and deodorizes. <0.15% alcohol content. Effective in 10 minutes or less against 100 microorganisms. Latex-free.

Canister of Extra Large Germicidal Surface Wipes. 65 wipes/can, 8" x 14" wipes.

| Number    | Quantity |
|-----------|----------|
| DI-3065   | 1/pkg    |
| DI-3065-6 | 6/case   |

Canister of Large Germicidal Surface Wipes. 160 wipes/can, 6" x 6.75" wipes.

| Number     | Quantity |
|------------|----------|
| DI-3160    | 1/pkg    |
| DI-3160-12 | 12/case  |

**WALL MOUNT SANI-BRACKET**

Adjustable Wall Mount Sani-Bracket for Large and Extra Large Sani-Cloth Wipes Canisters listed on this page (not compatible with CaviWipes canisters from page 105).

| Number  | Quantity |
|---------|----------|
| DI-1000 | 1/pkg    |



Sterile Processing Products – Disinfectants and Cleaners



Cavicide® disinfectant cleaner is deal for use on all non-porous surfaces. Non-caustic formula offers a 3 minute kill time for TB and a 2 minute kill time for HBV, HCV, HIV-1, MRSA, and VRE. Effective at removing mold and mildew. Safe on ultrasonic instruments. Fully biodegradable. Ready for use and requires no rinsing. Can be used in 3 different ways: an immersion solution, a surface disinfectant cleaner, or as an ultrasonic solution.

CAVIWIPES FLAT PACK  
CaviWipes Flat Pack. 45 wipes/pack, 7" x 9" wipes.

| Number     | Quantity |
|------------|----------|
| DI-4045    | each     |
| DI-4045-20 | 20/case  |

CAVIWIPES CANISTER  
Canister of CaviCide Extra Large Wipes. 66 wipes/can, 10" x 12" wipes.

| Number     | Quantity |
|------------|----------|
| DI-4066    | each     |
| DI-4066-12 | 12/case  |

Canister of CaviCide Large Wipes. 160 wipes/can, 6" x 6.75" wipes.

| Number     | Quantity |
|------------|----------|
| DI-4160    | each     |
| DI-4160-12 | 12/case  |

WALL BRACKET FOR CAVIWIPES CANISTER  
Wall Bracket for Large and Extra Large CaviWipes Canisters. 6" h x 5.5" w x 5.5" d (not compatible with SaniCloth canisters from page 103).

| Number  | Quantity |
|---------|----------|
| DI-4000 | each     |

CAVICIDE SPRAY BOTTLES  
CaviCide Surface Disinfectant 2 oz. Spray Bottle

| Number    | Quantity |
|-----------|----------|
| DI-5002   | each     |
| DI-5002-6 | 6/case   |

CaviCide Surface Disinfectant 8 oz. Spray Bottle

| Number     | Quantity |
|------------|----------|
| DI-5008    | each     |
| DI-5008-12 | 12/case  |

CaviCide Surface Disinfectant 24 oz. Spray Bottle

| Number     | Quantity |
|------------|----------|
| DI-5024    | each     |
| DI-5024-12 | 12/case  |

CaviCide Surface Disinfectant 1 Gallon Container

| Number    | Quantity |
|-----------|----------|
| DI-5128   | each     |
| DI-5128-4 | 4/case   |

CaviCide Surface Disinfectant with Spigot 2.5 Gallon Container

| Number    | Quantity |
|-----------|----------|
| DI-5256   | each     |
| DI-5256-2 | 2/case   |



# Sterile Processing Products – Enzymatic Detergents



A



B

**HEMAWAY ENZYMATIC FOAM®**  
 Apply on surgical instruments to pre-clean and prevent corrosion. Safe on metals, plastics and rubber. Non-hazardous, non-corrosive and biodegradable. Available by bottle or case.

| Number    | Description             | Details   | Quantity |
|-----------|-------------------------|-----------|----------|
| A IC-1116 | Hemaway Enzymatic Foam® | 16 ounces | 1 bottle |
| IC-1116-6 | Hemaway Enzymatic Foam® | 16 ounces | 6/case   |

**HEMAWAY ORIGINAL®**  
 Prolongs the life of instruments. Strong enough to remove the toughest debris, yet gentle enough to be used everyday. Safe on metals, plastics and rubber. Fresh, clean fragrance. Concentrated formula makes 128 gallons of solution. Non-hazardous, non-corrosive and biodegradable. Available by the bottle or by the case.

| Number    | Description       | Details  | Quantity |
|-----------|-------------------|----------|----------|
| B IC-1264 | Hemaway Original® | 1 gallon | 1 bottle |
| IC-1264-4 | Hemaway Original® | 1 gallon | 4/case   |

**HEMAWAY TRI-MAX®**  
 Greatly reduces the need for scrubbing and prolongs the life of instruments. Strong enough to remove the toughest debris, yet gentle enough to be used everyday. Safe on metals, plastics and rubber. Fresh, clean fragrance. Concentrated formula makes 128 gallons of solution. Non-hazardous, non-corrosive and biodegradable. Available by the bottle or by the case.

| Number    | Description      | Details  | Quantity |
|-----------|------------------|----------|----------|
| C IC-1364 | Hemaway Tri-Max® | 1 gallon | 1 bottle |
| IC-1364-4 | Hemaway Tri-Max® | 1 gallon | 4/case   |

**MPD® MULTI-PURPOSE DETERGENT**  
 Highly concentrated multi-purpose detergent formulated to clean surgical instruments gently and effectively. For manual or ultrasonic cleaning. Gentle enough to be used everyday. Safe on metals, plastics and rubber. Concentrated formula makes 128 gallons of solution. Non-hazardous, non-corrosive and biodegradable. Available by the bottle or by the case.

| Number    | Description             | Details  | Quantity |
|-----------|-------------------------|----------|----------|
| D IC-1464 | Multi-Purpose Detergent | 1 gallon | 1 bottle |
| IC-1464-4 | Multi-Purpose Detergent | 1 gallon | 4/case   |



C



D



## Sterile Processing Products – Efferzyme Cleaning Tablets, Wypalls™ and Polyester Wipers



### EFFERZYME – EFFERVESCENT CLEANING TABLETS

EfferZyme cleaning tablets are made from biodegradable, non-chlorine ingredients that contain fast acting protease enzymes. These neutral pH tablets enhance the removal of blood and other soils that may be resistant to chemical detergents, while still being delicate enough to not harm instruments or endoscopes. May be used as an ultrasonic cleaning solution or an instrument presoak. Removes the need to transport or store large jugs of liquid enzymes. Simply add two tablets per gallon of warm water.

| Number    | Quantity       |
|-----------|----------------|
| A IC-2152 | 52 tablets/box |

### KIMBERLY-CLARK WYPALL™ - L40 WIPERS

Kimberly-Clark Wypall™ - L40 Wipers are an outstanding general-purpose wiper featuring Kimberly-Clark's unique base sheet technology. Durable ¼ folded wipes easily absorb liquids yet are soft enough to use on delicate surfaces. Sold by the case.

| Number     | Description   | Quantity  |
|------------|---------------|---|
| B CR-21214 | 12.5" x 14.4" | 56 wipes/package, 18 packages/case, 1008 total/case |

### POLYESTER WIPER

Key Surgical Polyester Wiper is designed for cleanroom duties. Made from 100% continuous filament 85 denier medium weight polyester in a double-knit interlock, no run circular pattern to stop frays from starting. The wipes have ultrasonically sealed edges and are laundered and packaged in a controlled environment to reduce extractable particulates.

| Number     | Description | Quantity    |
|------------|-------------|-------------|
| C CR-20909 | 9" x 9"     | 150/package |

