



KEYSURGICAL®

KEY

8101 Wallace Road
Eden Prairie, MN 55344 USA

Tel. 800.541.7995 952.914.9789
Fax 952.914.9866

keysurgical.com info@keysurgical.com

KEY SURGICAL UTILITY SCISSORS AND OTHER GOOD STUFF

Key Surgical Good Stuff is where you'll find all the miscellaneous items that you just can't do without. We have all the good stuff needed to keep surgical and sterile processing departments running smoothly and efficiently.

Key Surgical Utility Scissors are designed with stainless steel blades that cut through just about everything. They're perfect for cutting bandages, casts, tubing and drapery, and they're equally at home in the O.R. and the E.R., in Central Sterile and in Recovery. Our selection of bright colors and two sizes just about guarantees that everyone will find their perfect pair.

CoFlex® NL was the first latex-free bandage to be truly cohesive. Our expanded line features sterile packs, kid prints, and camouflage CoFlex® NL patterns in variety of widths and colors. CoFlex® LF2 offers the next step forward in bandages with improved cohesiveness. CoFlex® AFD and CoFlex®-UMAFD intergrade absorbent bandages into CoFlex rolls for easier protection options. All CoFlex® options are easy to apply, don't constrict and tear clean easily - just three of the reasons why they're included in our Good Stuff section.

Key Surgical Good Stuff also includes an assortment of products that you can find elsewhere in this catalog. We've grouped things like Bundling Straps, Adhesive Remover, and our French Catheter Scale here for your convenience. We hope everyone can agree that this is all more Good Stuff from Key Surgical.

Utility Scissors

Bundling Straps

CoFlex® NL, CoFlex® LF2, CoFlex® AFD, and CoFlex®-UMAFD


Lint Removers

K-Wire Ruler and Pin Gauge

French Catheter Scale

Adhesive Remover

Brush Organizer





UTILITY SCISSORS

All-purpose shears with stainless steel blades and brightly colored handles. Excellent for cutting casts, bandages, tubing, and drapes. Designed for E.R., Central Sterile, Anesthesia, O.R., Orthopedics and general hospital use. Sold non-sterile. Autoclavable. Available in two sizes.



	Number	Color	Size
A	KS-065	magenta	6.5"
B	KS-066	neon orange	6.5"
C	KS-090	black	6.5"
D	KS-091	white	6.5"
E	KS-092	blue	6.5"
F	KS-093	green	6.5"
G	KS-093T	teal	6.5"
H	KS-094	red	6.5"
I	KS-095	yellow	6.5"
J	KS-096	Key Surgical blue	6.5"
K	KS-097	purple	6.5"
L	KS-098	lime green	6.5"
M	KS-099N	neon pink	6.5"
N	KS-100	black	7.5"
O	KS-101	white	7.5"
P	KS-102	blue	7.5"
Q	KS-103	green	7.5"
R	KS-103T	teal	7.5"
S	KS-104	red	7.5"
T	KS-105	yellow	7.5"
U	KS-106	Key Surgical blue	7.5"
V	KS-107	purple	7.5"
W	KS-108	lime green	7.5"
X	KS-109	neon pink	7.5"
Y	KS-165	magenta	7.5"
Z	KS-166	neon orange	7.5"





BUNDLING STRAPS

Hook and loop design allows for quick and easy bundling of electrical cords, IV lines and overhead microscope cords. Bright colors allow for quick identification of cords and lines. Designed for O.R., Recovery, E.R. and general hospital use. Latex-free. Sold non-sterile. 1" wide by 12" long. Thirty straps per package.

Number	Color
A MJ-402	blue
B MJ-403	green
C MJ-404	red
D MJ-405	yellow
E MJ-407	orange
F MJ-409	pink
MJ-400	assortment: five of each color

EMPLOYEE FEATURE

Jana

Job Title: **QA/RA Associate**

AKA: **Condimentaphobic**

What is your nickname? **Janban**. What is your guilty pleasure? **Online shopping and celebrity gossip**. What are you famous for around the office? **Being a picky eater and eating everything plain**. What is the last book you read? **Breaking Dawn**. What is your favorite place on earth? **Spider Lake in Hayward, Wisconsin on a nice summer day**. Name the one thing that you miss about being a kid. **Summer off and no worries**. When you were a child, what did you want to be when you grew up? **I wanted to be an orthopedic surgeon but I would never have survived med school**. What kind of car did you wish you had when you were a kid? **As a kid, a Ferrari because that is what my Barbie drove**. What were you doing the last time you really had a good laugh? **Reading textsfromlastnight.com**. What was the best vacation you ever had? **Going to London over New Year's Eve with my family**.



CoFlex® NL
 The first latex-free bandage to be truly cohesive, CoFlex NL is easy to apply, will not constrict and easily tears clean. Non-sterile CoFlex® NL assortment pack rolls are 5 yards in length and available in five widths.



Assorted Color Packs - Contain rolls of red, neon pink, neon green, light blue, blue, and purple.

Number	Size	Quantity
A CF-5100A	1" width	30/case
B CF-5200A	2" width	36/case
CF-5300A	3" width	24/case
CF-5400A	4" width	18/case
CF-5600A	6" width	12/case



Assorted Kid Color Packs - Contain rolls of orange smiles, yellow butterflies, pink hearts, purple flowers, blue cars, green dinosaurs.

Number	Size	Quantity
CF-5200KP	2" width	36/case
C CF-5300KP	3" width	24/case
CF-5400KP	4" width	18/case



Assorted Camouflage Color Packs - Contain rolls of pink, orange, green, olive green, blue, and purple camouflage.

Number	Size	Quantity
CF-5100CP	1" width	30/case
CF-5200CP	2" width	36/case
CF-5300CP	3" width	24/case
D CF-5400CP	4" width	18/case
CF-5600CP	6" width	12/case

PRODUCTS SHOWN 50% ACTUAL SIZE



PRODUCTS SHOWN ACTUAL SIZE

CoFlex® NL
 The first latex-free bandage to be truly cohesive, CoFlex® NL is easy to apply, will not constrict and easily tears clean. Non-sterile CoFlex® NL rolls are five yards in length and available in five widths.

Single Color Packs - Tan (Also available in Red, Blue, Green, Purple and other colors. Please call for more information.)

Number	Size	Quantity
A CF-5100T	1" width	30/case
B CF-5200T	2" width	36/case
C CF-5300T	3" width	24/case
D CF-5400T	4" width	18/case
E CF-5600T	6" width	12/case

CoFlex® NL Sterile Packs
 The first latex-free bandage to be truly cohesive, CoFlex® NL is easy to apply, will not constrict and easily tears clean. Sterile CoFlex® NL rolls are five yards in length and available in three widths and one color (tan). Sold sterile.

Number	Size	Quantity
F CF-5300S	3" width	24/case
CF-5400S	4" width	18/case
CF-5600S	6" width	12/case





CoFlex® LF2
 The next generation of CoFlex® latex-free cohesive bandage. CoFlex® LF2 is easy to apply, will not constrict, easily tears clean, and has improved cohesiveness. Non-sterile CoFlex® LF2 rolls are five yards in length and available in five widths.

Assorted Color Packs - Contain rolls of red, neon pink, neon green, light blue, blue, and purple.

Number	Size	Quantity
A CF-9100A	1" width	30/case
B CF-9200A	2" width	36/case
C CF-9300A	3" width	24/case
CF-9400A	4" width	18/case
CF-9600A	6" width	12/case

Single Color Packs - Tan (Also available in red, blue, green, purple and other colors. Please call for more information.)

Number	Size	Quantity
D CF-9100T	1" width	30/case
E CF-9200T	2" width	36/case
F CF-9300T	3" width	24/case
CF-9400T	4" width	18/case
CF-9600T	6" width	12/case

PRODUCTS SHOWN 50% ACTUAL SIZE



CoFlex® AFD
 A latex-free cohesive bandage with Absorbent Foam Dressing layer. All-in-one emergency cohesive bandage system. Foam pad absorbs and retains fluid. Waterproof bandage layer prevents leakage and contamination. 4" width by 2.5 yards, available in sterile and non-sterile.

Number	Details	Quantity
A CF-9740T	tan, non-sterile	8/case
B CF-9740S	tan, sterile	8/case

CoFlex® + UMAFD
 A latex-free cohesive bandage with Universal Medic Absorbent Foam layer. All-in-one cohesive trauma bandage system. Absorbent foam pad contains waterproof film layer to absorb fluid and prevent leaks. Vacuum packaged. 4" width by 3 yards. Non-sterile.

Number	Details	Quantity
C CF-9438UM	black	18/case
CF-9430UM	black	100/case

EMPLOYEE FEATURE

Kasey
 Job Title: **QA/RA Associate**
 AKA: **Chicken Scratch**

What is your sign? **Libra.** What is your favorite place on earth? **The lake on a calm day.** What is usually your first thought when you wake up? **Coffee!!** What don't you like about living in Minnesota? **That I've been here my entire life.** Name three exotic countries that you would like to visit. **Thailand, New Zealand, Chile.** If you could have any super power, what would it be? **Time travel, it could definitely be useful.** What is your favorite flavor of jelly bean? **I don't like jelly beans – weird, huh?** What song always makes you cringe and turn off the music? **Anything country.** What is your favorite office machine? **The phones – I love calling people who are sitting right next to me.** What is the last magazine you read? **Runner's World.** What are you famous for around the office? **Illegible handwriting.** What is your favorite thing to do on a Sunday afternoon? **Nothing.**

Good Stuff – Lint Removers, K-Wire Ruler and Pin Gauge, French Catheter Scale, Adhesive Remover and Brush Organizer



LINT REMOVER HANDLE AND LINT REMOVER REFILLS

Designed for heavy-duty hospital use.

	Number	Dimensions	Quantity
A	MJ-001	9" overall length	1 handle
B	MJ-012	60 sheets per roll	12 rolls

K-WIRE RULER AND PIN GAUGE

Practical tool measures inches, millimeters and 16 sizes of wires and pins. Made from 420-series stainless steel with .025" thickness for durability.

	Number	Details	Quantity
C	KI-39-710	6" ruler	1

FRENCH CATHETER SCALE

Convenient two-sided gauge allows the right choice to clean each instrument.

	Number	Description	Quantity
D	MJ-370	plastic card with die-cut holes	4

ADHESIVE REMOVER

Specially formulated to be gentle on surgical instruments while effectively removing tape and adhesive residue. Restores surface to original condition. Fresh citrus fragrance. Environmentally safe. No hazardous solvents or toxic fumes. Available in two sizes.

	Number	Details	Quantity
E	MJ-7004	4 ounces	box of six bottles
F	MJ-7032	32 ounces each	one bottle

BRUSH ORGANIZER

Stainless steel organizer mounts easily on wall and can be customized to fit individual needs. 12 permanent pegs. 2 adjustable bars with 4 pegs per bar. 16" wide by 32" long.

	Number	Description
G	BR-7000	wall-mounted organizer

BRUSH ORGANIZER WITH BRUSH ASSORTMENT (BR-4199)

Assortment of brushes includes one package each of the following brushes:

BR-18-062	BR-18-087	BR-18-105
BR-18-125	BR-18-250	FT-18-500
N-2100	N-2121	S-2100

Number	Description
BR-7100	wall-mounted organizer with brush assortment

ADJUSTABLE BAR FOR BRUSH ORGANIZER

Extra adjustable bar with 4 pegs.

	Number	Description
H	BR-7050	adjustable bar with four pegs